


EP Tank Biosense Webster

PACKING GUIDE



LIST OF COMPONENTS:

- A** EP Tank table (#1)
- B** EP Tank (#1)
- C** Power Supply (#1)
- D** Power Supply Cables (#2): Cable EP Tank-Power Supply and Cable Power Supply- AC outlet
- E** Table legs (#5)
- F** Catheters support (#1)
- G** Mini table (#1)
- H** Endoscope (#1)
- I** External pump (to use to empty the EP Tank) (#1) - *not supplied after version 2022** -
- L** Tube for the external pump (#1)
- M** Screwdriver (#1)
- N** Towel to clean the EP Tank after the use (#1)
- O** Pins (#10)

**Please refer to the "Cleaning Guide" to know how to clean the EP-Tank*